

MINI MASSAGER

Relaxation ,Anywhere, Anytime

Included in this standard package:

- Device x 1pc
- Large pad x 2pcs
- Small pad x 2pcs
- Electric wire x 2pcs
- USB cable & AC adapter x 1set
- Acupuncture point chart x 1pc
- Instruction manual x 1pc
- Collection board for pad and wire x 1pc


Combo package also available:

- XL electrode pad
- Shoes for feet stimulation

Features

- 2 output with 6 modes
- 20 level intensity
- 20 minutes default timer, can adjustable from 10-60 minutes
- Rechargeable lithium battery
- Battery indicator to know when need to recharge battery
- AC adapter 110-240v, USB Charger 5v

Illustrations of the device and accessories: (Model:AS1088)



- | | |
|---------------------------------|--------------------------------------|
| 1. Two output connector | 8. Time adjustment Battery indicator |
| 2. On/ off switch | 9. Pads- sticking board |
| 3. USB port | 10. Massage pads |
| 4. LCD screen | 11. Pad wires |
| 5. Mode selector | 12. USB cable |
| 6. Intensity adjustment | 13. Charger |
| 7. Pause switch | |

Modes

There are six modes to choose from, there are:

1. Pushing
2. Acupuncture
3. Tapping
4. Scraping
5. Cupping
6. Random (This mode cycles all 5 modes)

Operating instructions

1. Connect the pads to the wire by snapping them on, then remove the protective film.

Note: adding a few drops of water on the pads will increase their longevity.

2. Place pads where the massage sensations are desired. And plug into one of the ports on the device.

For additional relief use the second port and repeat procedure to put on another two pads.

Make sure both pads on the skin and not overlapping. Make sure your skin is free from any dirt, oil or lotion.

3. Turn the device on by sliding the on/ off switch on the bottom of the unit. The unit will be now activated.

4. Press the mode button to try different mode. Remember the strength will default to the lowest setting by each mode selection.

Note: Center button can also be used as a lock key.

Press it to lock the screen, again to unlock the screen.

5. Once you select a mode, gradually increase the intensity by pressing the + button; decrease it by pressing - button.

6. Recommended application is 10-60 minutes. Increase 10 minutes by each setting. The time remaining is automatically counted down. Then the unit will automatically turn off after time

finish.

Precautions

- Do not apply the massage pads near the heart, on the head, above the neck, in the pubic region, over Scarred areas, on the throat or over the mouth.
- Avoid touching the pads when the unit is on.
- Do not use this massager if you a pacemaker or other life support equipment such as an artificial lung or respirator.
- should not be used while driving, operating machinery, or during any activity in which involuntary muscle contractions may put the user at undue risk or injury
- Do not use if you have following medical conditions:

Acute disease, trauma or surgical procedure in past six month

Cardiac dysfunction

Epilepsy

High blood pressure

High fever

cancer

Infectious disease

Malignant tumors

Recent surgery

Skin problems

Undiagnosed pain syndromes

Abdominal or inguinal hernia

Painful and /or atrophied muscles

Limited range of motion in skeletal joints

- Keep out of the reach of children.
- the unit is not intended for the application of any medical condition or disease nor is it intended for physiotherapy or muscle rehabilitation. It is contraindicated for use on any muscle that is injured or diseased

Maintenance and care

1. When you are finished using the device, turn it off and take off the pads. Store the unit in a cool, dry place until next use.
2. After usage, please cover the pads with the protective film and put them back in their zip lock bags. Applying a few drops water before and after usage will increase the longevity of the pads.
3. Charge the unit at least every 3 month in event it is not used regularly. You can also use the unit with the wall charger or directly from your computer using the USB cable.

Recommend massage points



Troubleshooting

Problem	Cause	Solution
One pad feels stronger than the other.	This is normal. Different areas of your body will react differently.	Nothing needs to be done. Make sure the pads are moist and are making good contact.
During the massage, the skin feels a painful burning sensation or the stimulation becomes weakened.	Adhesive gel pads are not adhering firmly to the skin. Or the gel pads are too dry.	Apply a few drops of water to adhesive surface of each pad and make sure the pads are pressed firmly to the skin during application.
Unit is on, but no	Both pads are not	Make sure both pads are

sensation is felt in the pads.	firmly on the skin. Or there is a loose connection. Or the intensity level needs to be increased.	firmly pressed to the skin. Make sure all connections are secure from the unit to the wires and the wires to the pads. Increase the intensity level.
Adhesive gel pads do not stick to skin even after cleaning and moistening the gel pad.	Adhesive gel pads need to be replaced.	Replace the set of pads. Pads can be purchased from local distributor.
Unit does not turn on or the screen is dim.	Battery is low.	Recharge the battery.
Does not provide adequate pain relief.	Not using it long enough. Not lying down then massaging neck.	Use your massager 20-30 minutes at a time, 3-6 times per day. Lie on your back for neck support when you need relief for headaches or neck pain.